



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KAPCB Antibody [#C10453]

Catalog Number: C10453

Concentration: 1mg/ml

Swiss-Prot No.: P22694

Other Names: EC 2.7.11.11, KAPB, KAPCB, PKA C-beta, PKA-beta, cAMP-dependent protein kinase, beta-catalytic subunit

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

KAPCB.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAPCB antibody detects endogenous levels of total KAPCB protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Beebe S.J., Mol. Endocrinol. 4:465-475(1990).

Gregory S.G., Nature 441:315-321(2006).

Wu K.-J., Oncogene 21:7872-7882(2002).

