



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SFRS4 Antibody [#C10906]

Catalog Number: C10906

Concentration: 1mg/ml

Swiss-Prot No.: Q08170

Other Names: Splicing factor, arginine/serine-rich 4, Pre-mRNA-splicing factor SRP75, SRP001LB, SFRS4, SRP75

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

SFRS4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SFRS4 antibody detects endogenous levels of total SFRS4 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Zahler A.M., Mol. Cell. Biol. 13:4023-4028(1993).

Gregory S.G., Nature 441:315-321(2006).

Zahler A.M., Genes Dev. 6:837-847(1992).

