



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MEKKK 1 Antibody [#C11444]

Catalog Number: C11444

Concentration: 1mg/ml

Swiss-Prot No.: Q92918

Other Names: Mitogen-activated protein kinase kinase kinase 1, EC 2.7.11.1, MAPK/ERK kinase kinase 1, MEK kinase kinase 1, MEKKK 1, Hematopoietic progenitor kinase, MAP4K1, HPK1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MEKKK 1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MEKKK 1 antibody detects endogenous levels of total MEKKK 1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Hu M.C.-T., Genes Dev. 10:2251-2264(1996).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).

Wissing J., Mol. Cell. Proteomics 6:537-547(2007).

