



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

TF3B Antibody [#C11452]

Catalog Number: C11452

Concentration: 1mg/ml

Swiss-Prot No.: Q92994

Other Names: TFIIB90; hTFIIB90; B-related factor 1; BRF-1; hBRF; TATA box-binding protein-associated factor, RNA polymerase III, subunit 2; TAF3B2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TF3B.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TF3B antibody detects endogenous levels of total TF3B protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
 IF: 1:100~1:500 ELISA: 1:1000

References:

Wang Z., Proc. Natl. Acad. Sci. U.S.A. 92:7026-7030(1995).

Mital R., Mol. Cell. Biol. 16:7031-7042(1996).

McCulloch V., EMBO J. 19:4134-4143(2000).





Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TF3B antibody.



Peptide - - +

Western blot analysis of extracts from HeLa cells and Jurkat cells, using TF3B antibody.



Peptide

-

+

Immunofluorescence analysis of MCF-7 cells, using TF3B antibody.