



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CIB3 Antibody [#C11530]

Catalog Number: C11530

Concentration: 1mg/ml

Swiss-Prot No.: Q96Q77

Other Names: Calcium and integrin-binding family member 3, Kinase-interacting protein 3, KIP 3, CIB3, KIP3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human CIB3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CIB3 antibody detects endogenous levels of total CIB3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Hayashi A., Submitted (NOV-2000) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

