



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ALPK1 Antibody [#C11534]

Catalog Number: C11534

Concentration: 1mg/ml

Swiss-Prot No.: Q96QP1

Other Names: ALPK1; EC 2.7.11.-; LAK; Lymphocyte alpha-kinase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ALPK1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ALPK1 antibody detects endogenous levels of total ALPK1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ryazanov A.G., Curr. Biol. 9:R43-R45(1999).

Nagase T., DNA Res. 7:143-150(2000).

Yamada S., Oncogene 23:5901-5911(2004).

