



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARSK Antibody [#C14573]

Catalog Number: C14573

Concentration: 1mg/ml

Swiss-Prot No.: Q6UWY0

Other Names: Arylsulfatase K; EC 3.1.6.-; ASK; Telethon sulfatase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ARSK.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARSK antibody detects endogenous levels of total ARSK protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 94%, Positives = 100%);
Rat (Identities = 88%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Sardiello M., Hum. Mol. Genet. 14:3203-3217(2005).

Clark H.F., Genome Res. 13:2265-2270(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using ARSK antibody.