



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

Collagen XX α 1 Antibody [#C12225]

Catalog Number: C12225

Concentration: 1mg/ml

Swiss-Prot No.: Q9P218

Other Names: Collagen alpha-1(XX) chain [Precursor], COKA1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human Collagen XX α 1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Collagen XX α 1 antibody detects endogenous levels of total Collagen XX α 1 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
 IF: 1:100~1:500 ELISA: 1:20000

References:

Deloukas P., Nature 414:865-871(2001).


Nagase T., DNA Res. 7:143-150(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using Collagen XX α 1 antibody.



Peptide - +

Western blot analysis of extracts from HUVEC cells, using Collagen XX α 1 antibody.



Peptide - +

Immunofluorescence analysis of NIH/3T3 cells, using Collagen XX α 1 antibody.

