



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Cytochrome P450 3A43 Antibody [#C12274]

**Catalog Number:** C12274

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9HB55

**Other Names:** Cytochrome P450 3A4, EC 1.14.14.1, CYP3A43

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 3A43.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Cytochrome P450 3A43 antibody detects endogenous levels of total Cytochrome P450 3A43 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Domanski T.L., Mol. Pharmacol. 59:386-392(2001).

Westlind A., Biochem. Biophys. Res. Commun. 281:1349-1355(2001).

Gellner K., Pharmacogenetics 11:111-121(2001).

