



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IL17RAntibody [#C12431]

Catalog Number: C12431

Concentration: 1mg/ml

Swiss-Prot No.: Q96F46

Other Names: IL-17 receptor; IL-17RA; IL17RA;
interleukin-17 receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IL17RA.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IL17RAntibody detects endogenous levels of total IL17RAprotein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 92%, Positives = 100%)

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:60000

References:

Dunham I., Genome Biol. 5:R84.1-R84.11(2004).

Venter J.C., Submitted (JUL-2005).

Henzel W.J., Protein Sci. 13:2819-2824(2004).

