



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NF1 Antibody [#C13094]

Catalog Number: C13094

Concentration: 1mg/ml

Swiss-Prot No.: P21359

Other Names: Neurofibromatosis-related protein NF-1;
Neurofibromin

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human NF1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: NF1 antibody detects endogenous levels of
total NF1 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Bernards A., DNA Cell Biol. 11:727-734(1992).

Wallace M.R., Science 249:181-186(1990).

Wallace M.R., Science 250:1749-1749(1990).



Western blot analysis of extracts from HepG2 cells, using
NF1 antibody.



Immunofluorescence analysis of HepG2 cells, using NF1
antibody.