



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PHLA1 Antibody [#C13101]

Catalog Number: C13101

Concentration: 1mg/ml

Swiss-Prot No.: Q8WV24

Other Names: Pleckstrin homology-like domain family A member 1; T-cell death-associated gene 51 protein; Apoptosis-associated nuclear protein; Proline- and histidine-rich protein; Proline- and glutamine-rich protein; PQ-rich protein; PHLDA1; PHRIP, TDAG51

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human PHLA1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PHLA1 antibody detects endogenous levels of total PHLA1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Kuske M.D., Cytogenet. Cell Genet. 89:1-1(2000).

Neef R., Cancer Res. 62:5920-5929(2002).

Hinz T., Cell. Signal. 13:345-352(2001).

