



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

RPL17 Antibody [#C14153]

Catalog Number: C14153

Concentration: 1mg/ml

Swiss-Prot No.: P18621

Other Names: 60S ribosomal L17; gene encoding putative NFkB activating protein; L23; ribosomal protein L17; RL17; rpl23

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RPL17.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPL17 antibody detects endogenous levels of total RPL17 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
 ELISA: 1:10000

References:

Mager D.L., Nucleic Acids Res. 18:5301-5301(1990).


Batra S.K., Cell Growth Differ. 2:385-390(1991).

Yoshihama M., Genome Res. 12:379-390(2002).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue using RPL17 antibody.



Western blot analysis of extracts from HeLa cells and COS cells, using RPL17 antibody.