



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CG028 Antibody [#G223]

Catalog Number: G223

Concentration: 1mg/ml

Swiss-Prot No.: O95766

Other Names: UPF0550 protein C7orf28; C7orf28A; CGI-43; C7orf28B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CG028.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CG028 antibody detects endogenous levels of total CG028 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Lai C.-H., Genome Res. 10:703-713(2000).

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 424:157-164(2003).

