



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Maf Antibody [#C10174]

Catalog Number: C10174

Concentration: 1mg/ml

Swiss-Prot No.: O75444

Other Names: Transcription factor Maf, Proto-oncogene c-maf, MAF

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human Maf.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Maf antibody detects endogenous levels of total Maf protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Chesi M., Blood 91:4457-4463(1998).

Jamieson R.V., Hum. Mol. Genet. 11:33-42(2002).

Vanita V., Am. J. Med. Genet. A 140:558-566(2006).

