



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GALR3 Antibody [#G269]

Catalog Number: G269

Concentration: 1mg/ml

Swiss-Prot No.: O60755

Other Names: Galanin receptor type 3; GAL3-R; GALR3; GALR3; GALNR3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GALR3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GALR3 antibody detects endogenous levels of total GALR3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Smith K.E., J. Biol. Chem. 273:23321-23326(1998).

Bennett M.M., Submitted (MAY-1998) to the EMBL/GenBank/DDBJ databases.

Kolakowski L.F. Jr., J. Neurochem. 71:2239-2251(1998).

