



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TACD2 Antibody [#C10320]

Catalog Number: C10320

Concentration: 1mg/ml

Swiss-Prot No.: P09758

Other Names: EGP-1, GA733-1, M1S1, TACD2, TROP2, cell surface glycoprotein Trop-2, pancreatic carcinoma marker protein GA733-1, tumor-associated calcium signal transducer 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

TACD2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TACD2 antibody detects endogenous levels of total TACD2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Linnenbach A.J., Proc. Natl. Acad. Sci. U.S.A. 86:27-31(1989).

Fornaro M., Int. J. Cancer 62:610-618(1995).

Tsujikawa M., Nat. Genet. 21:420-423(1999).

