



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MYF5 Antibody [#C10359]

Catalog Number: C10359

Concentration: 1mg/ml

Swiss-Prot No.: P13349

Other Names: MYF5, myogenic factor 5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MYF5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MYF5 antibody detects endogenous levels of total MYF5 protein.


Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Braun T., EMBO J. 8:701-709(1989).

Thompson D.B.; Submitted (MAY-1995) to the EMBL/GenBank/DDBJ databases.



For Research Use Only