



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## hnRNP A2/B1 Antibody [#C10452]

**Catalog Number:** C10452

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P22626

**Other Names:** HNRPA2B1, Heterogeneous nuclear ribonucleoproteins A2/B1, ROA2, hnRNP A2 / hnRNP B1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human hnRNP A2/B1.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** hnRNP A2/B1 antibody detects endogenous levels of total hnRNP A2/B1 protein.

**Reactivity:** Human, Mouse


**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:20000

**References:**


Burd C.G., Proc. Natl. Acad. Sci. U.S.A. 86:9788-9792(1989).

Kozu T., Genomics 25:365-371(1995).

Biamonti G., Nucleic Acids Res. 22:1996-2002(1994).



Western blot analysis of extracts from HepG2 cells, using hnRNP A2/B1 antibody.



Immunofluorescence analysis of HeLa cells, using hnRNP A2/B1 antibody.