



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IκB-α Antibody [#C10481]

Catalog Number: C10481

Concentration: 1mg/ml

Swiss-Prot No.: P25963

Other Names: I-kappa-B-alpha; IkappaBalpha; IKBA; MAD3; Major histocompatibility complex enhancer-binding protein MAD3; NF-kappaB inhibitor alpha; NFKBI; NFKBIA; RL/IF-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human IκB-α.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IκB-α antibody detects endogenous levels of total IκB-α protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Haskill S., Cell 65:1281-1289(1991).

Jungnickel B., J. Exp. Med. 191:395-402(2000).

Liu B., Submitted (APR-2001) to the EMBL/GenBank/DDBJ databases.



For Research Use Only