



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR171 Antibody [#G310]

Catalog Number: G310

Concentration: 1mg/ml

Swiss-Prot No.: O14626

Other Names: Probable G-protein coupled receptor 171;
G-protein coupled receptor H963; GPR171; H963

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

GPR171.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: GPR171 antibody detects endogenous levels
of total GPR171 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Jacobs K.A., Gene 198:289-296(1997).

Joensuu T., Am. J. Hum. Genet. 69:673-684(2001).

