



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## EMR3 Antibody [#G092]

**Catalog Number:** G092

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BY15

**Other Names:** EGF-like module-containing mucin-like hormone receptor-like 3 [Precursor]; EGF-like module-containing mucin-like hormone receptor-like 3 [Precursor]; EGF-like module-containing mucin-like hormone receptor-like 3 subunit alpha; EGF-like module-containing mucin-like hormone receptor-like 3 subunit beta; EMR3; UNQ683/PRO1562

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human EMR3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** EMR3 antibody detects endogenous levels of total EMR3 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Stacey M., J. Biol. Chem. 276:18863-18870(2001).

Clark H.F., Genome Res. 13:2265-2270(2003).

Olsen J.V., Cell 127:635-648(2006).

