



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FZD10 Antibody [#G106]

Catalog Number: G106

Concentration: 1mg/ml

Swiss-Prot No.: Q9ULW2

Other Names: Frizzled-10 [Precursor]; Fz-10; hFz10; FzE7; CD350 antigen; FZD10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FZD10.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FZD10 antibody detects endogenous levels of total FZD10 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Koike J., Biochem. Biophys. Res. Commun. 262:39-43(1999).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Tanaka S., Proc. Natl. Acad. Sci. U.S.A. 95:10164-10169(1998).

