



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PITX1 Antibody [#C10791]

Catalog Number: C10791

Concentration: 1mg/ml

Swiss-Prot No.: P78337

Other Names: Pituitary homeobox 1, Hindlimb-expressed homeobox protein backfoot, PITX1, BFT, PTX1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PITX1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PITX1 antibody detects endogenous levels of total PITX1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Shang J., Genomics 40:108-113(1997).

Schmutz J., Nature 431:268-274(2004).

Olsen J.V., Cell 127:635-648(2006).

