



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ST5 Antibody [#C10797]

**Catalog Number:** C10797

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P78524

**Other Names:** Suppression of tumorigenicity 5, HeLa tumor suppression 1, DENN domain-containing protein 2B, ST5, DENND2B, HTS1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ST5.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ST5 antibody detects endogenous levels of total ST5 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Lichy J.H., Nucleic Acids Res. 24:4700-4708(1996).

Amid C., Cytogenet. Cell Genet. 93:284-290(2001).

Lichy J.H., Cell Growth Differ. 3:541-548(1992).

