



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

CRBP III Antibody [#C10805]

Catalog Number: C10805

Concentration: 1mg/ml

Swiss-Prot No.: P82980

Other Names: Retinol-binding protein III, cellular, CRBP-III, HRBPiso, RBP5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CRBP III.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CRBP III antibody detects endogenous levels of total CRBP III protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
 ELISA: 1:20000

References:

Durkin J.T., Biochemistry 43:16348-16355(2004).


The MGC Project Team, Genome Res. 14:2121-2127(2004).

McNee J.J., Submitted (APR-2000) to the EMBL/GenBank/DDBJ databases.



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using CRBP III antibody.



Peptide - +

Immunofluorescence analysis of LOVO cells, using CRBP III antibody.