



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

TNAP2 Antibody [#C10867]

Catalog Number: C10867

Concentration: 1mg/ml

Swiss-Prot No.: Q03169

Other Names: Tumor necrosis factor, alpha-induced protein 2, Primary response gene B94 protein, TNFAIP2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TNAP2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TNAP2 antibody detects endogenous levels of total TNAP2 protein.


Reactivity: Human (Identities = 100%, Positives = 100%);
 Mouse (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 IF: 1:100~1:500
 ELISA: 1:5000

References:


Sarma V., J. Immunol. 148:3302-3312(1992).

Rieder M.J., Submitted (APR-2003) to the EMBL/GenBank/DDBJ databases.



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue, using TNAP2 antibody.



Peptide - +

Immunofluorescence analysis of HuvEc cells, using TNAP2 antibody.