



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FRK Antibody [#C10609]

Catalog Number: C10609

Concentration: 1mg/ml

Swiss-Prot No.: P42685

Other Names: EC 2.7.1.112, EC 2.7.10.2, Nuclear tyrosine protein kinase RAK, Tyrosine-protein kinase FRK, fyn-related kinase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FRK.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FRK antibody detects endogenous levels of total FRK protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
IF: 1:100~1:500 ELISA: 1:1000

References:

Lee J., Gene 138:247-251(1994).


Cance W.G., Cell Growth Differ. 5:1347-1355(1994).

Zhang Y., Mol. Cell. Proteomics 4:1240-1250(2005).




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using FRK antibody.



Peptide - +

Western blot analysis of extracts from Jurkat cells, using FRK antibody.



Peptide - +

Immunofluorescence analysis of HUVEC cells, using FRK antibody.