



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Flt3 ligand Antibody [#C10654]

**Catalog Number:** C10654

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P49771

**Other Names:** SL cytokine [Precursor], Fms-related tyrosine kinase 3 ligand, Flt3 ligand, Flt3L, FLT3LG

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Flt3 ligand.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Flt3 ligand antibody detects endogenous levels of total Flt3 ligand protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Hannum C., Nature 368:643-648(1994).

Lyman S.D., Blood 83:2795-2801(1994).

Savvides S.N., Nat. Struct. Biol. 7:486-491(2000).

