



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BRSK1 Antibody [#C11392]

Catalog Number: C11392

Concentration: 1mg/ml

Swiss-Prot No.: Q8TDC3

Other Names: EC 2.7.11.1, SAD1 kinase,
Serine/threonine kinase SAD-B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BRSK1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BRSK1 antibody detects endogenous levels of total BRSK1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Lizcano J.M., EMBO J. 23:833-843(2004).

Lu R., J. Biol. Chem. 279:31164-31170(2004).

Nagase T., DNA Res. 8:85-95(2001).

