



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AP2C Antibody [#C11429]

Catalog Number: C11429

Concentration: 1mg/ml

Swiss-Prot No.: Q92754

Other Names: AP2-gamma, Activating enhancer-binding protein 2 gamma, Transcription factor ERF-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AP2C.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AP2C antibody detects endogenous levels of total AP2C protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Williamson J.A., Genomics 35:262-264(1996).

McPherson L.A., Proc. Natl. Acad. Sci. U.S.A. 94:4342-4347(1997).

Deloukas P., Nature 414:865-871(2001).



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using AP2C antibody.