



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TAF15 Antibody [#C11436]

Catalog Number: C11436

Concentration: 1mg/ml

Swiss-Prot No.: Q92804

Other Names: RNA-binding protein 56, TAFII68, TAF(II)68

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TAF15.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TAF15 antibody detects endogenous levels of total TAF15 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:1000

References:

Morohoshi F., Genomics 38:51-57(1996).


Bertolotti A., EMBO J. 15:5022-5031(1996).

Morohoshi F., Gene 221:191-198(1998).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TAF15 antibody.



Peptide - +

Western blot analysis of extracts from HepG2 cells, using TAF15 antibody.