



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MIC1 Antibody [#C11486]

Catalog Number: C11486

Concentration: 1mg/ml

Swiss-Prot No.: Q96DM3

Other Names: Uncharacterized protein C18orf8; Colon cancer-associated protein Mic1; Mic-1; C18orf8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human MIC1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: MIC1 antibody detects endogenous levels of total MIC1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Linnebacher M., Submitted (APR-1999) to the EMBL/GenBank/DDBJ databases.

