



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OLFML2A Antibody [#G484]

Catalog Number: G484

Concentration: 1mg/ml

Swiss-Prot No.: Q68BL7

Other Names: olfactomedin-like; olfactomedin-like 2A; OLFML2A; photomedin-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OLFML2A.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OLFML2A antibody detects endogenous levels of total OLFML2A protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Clark H.F., Genome Res. 13:2265-2270(2003).

Humphray S.J., Nature 429:369-374(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from RAW264.7 cells, using OLFML2A antibody.



Immunofluorescence analysis of A549 cells, using OLFML2A antibody.