



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2M7 Antibody [#G563]

Catalog Number: G563

Concentration: 1mg/ml

Swiss-Prot No.: Q8NG81

Other Names: Olfactory receptor 2M7; Olfactory receptor OR1-58

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR2M7.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2M7 antibody detects endogenous levels of total OR2M7 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Gregory S.G., Nature 441:315-321(2006).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).



Western blot analysis of extracts from HepG2 cells, Jurkat cells, MCF-7 cells, COLO cells and HUVEC cells, using OR2M7 antibody.



Immunofluorescence analysis of MCF-7 cells, using OR2M7 antibody.