



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## OR8G1 Antibody [#G686]

**Catalog Number:** G686

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q15617

**Other Names:** HSTPCR25; olfactory receptor 8G1;  
olfactory receptor OR11-281; olfactory receptor TPCR25;  
olfactory receptor, family 8, subfamily G, member 1;  
olfactory receptor, family 8, subfamily G, member 1  
pseudogene; OR8G1P; TPCR25

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against  
synthesized peptide derived from C-terminal of human

OR8G1.

**Purification:** The antibody was affinity-purified from  
rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** OR8G1 antibody detects endogenous levels  
of total OR8G1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Suwa M., Submitted (JUL-2001) to the  
EMBL/GenBank/DDBJ databases.

Vanderhaeghen P., Genomics 39:239-246(1997).

Fuchs T., Genomics 80:295-302(2002).

