



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CATL2 (Cleaved-Leu114) Antibody [#L0185]

Catalog Number: L0185

Concentration: 1mg/ml

Swiss-Prot No.: O60911

Other Names: Cathepsin L2 [Precursor], EC 3.4.22.43, Cathepsin V, Cathepsin U, CTSL2, CATL2, CTSU, CTSV, UNQ268/PRO305

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human CATL2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CATL2 (Cleaved-Leu114) antibody detects endogenous levels of fragment of activated CATL2 resulting from cleavage adjacent to Leu114.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Santamaria I., Cancer Res. 58:1624-1630(1998).

Adachi W., Invest. Ophthalmol. Vis. Sci. 39:1789-1796(1998).

Broemme D., Biochemistry 38:2377-2385(1999).



For Research Use Only