



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

COX17 Antibody [#C12231]

Catalog Number: C12231

Concentration: 1mg/ml

Swiss-Prot No.: Q14061

Other Names: Cytochrome c oxidase copper chaperone, COX17

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human COX17.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: COX17 antibody detects endogenous levels of total COX17 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
 ELISA: 1:5000

References:

Amaravadi R., Hum. Genet. 99:329-333(1997).


Punter F.A., Hum. Genet. 107:69-74(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using COX17 antibody.



Western blot analysis of extracts from HeLa cells, HepG2 cells, COLO cells and 293 cells, using COX17 antibody.