



Web: www.anbobio.com
www.dijibio.com
 Order: order@dijibio.com
 Tel: (86) 519 8805 0026

CIDEB Antibody [#C13034]

Catalog Number: C13034

Concentration: 1mg/ml

Swiss-Prot No.: Q9UHD4

Other Names: cell death activator CIDE-B; cell death-inducing DFFA-like effector B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CIDEB.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CIDEB antibody detects endogenous levels of total CIDEB protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
 IF: 1:100~1:500 ELISA: 1:40000

References:

Lugovskoy A.A., Cell 99:747-755(1999).


Inohara N., EMBO J. 17:2526-2533(1998).

Ota T., Nat. Genet. 36:40-45(2004).




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human cervix carcinoma tissue using CIDEB antibody.



Peptide - - +

Western blot analysis of extracts from RAW264.7 cells and HT-29 cells, using CIDEB antibody.



Peptide - +

Immunofluorescence analysis of COS-7 cells, using CIDEB antibody.